

SSAFlab Filtration Kit for ASTM D5452, IP 423 19722-0

ASTM D5452; IP 423

SSAFlab Filtration Kit for the determination of particulate contamination in aviation fuels.

Comprises 0.45 micron Filter, matched weight pair of 0.8 micron Filters, stainless steel Filter Holder with stainless steel screen, Filter Forceps, 25mm Solvent Filtering Dispenser, Vacuum Filtering Flask, and Petrislide.



Filtration System for ASTM D6217, IP 415 19724-0

ASTM D6217; IP 415

Filtration system for testing middle distillates in accordance with ASTM D6217 and IP 415.

Comprises 0.8 micron Filter, stainless steel Filter Holder with stainless steel screen, Filter Forceps, 25mm Solvent Filtering Dispenser, 2 x Vacuum Filtering Flasks 1 litre cylinder, Watch Glass and Tubing.

